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
RUBBER-DAM set

RUBBER DAM Instruments

BEFORE USE, READ THE FOLLOWING INSTRUCTION CAREFULLY

FOR DENTAL USE ONLY

NON-STERILE - STERILIZE DIRECTLY BEFORE USE

SET INCLUDED:

1. Profiled rubber dam frames available in 3 sizes.
2. Rubber dam forceps for clamps placement with the opening degree resistance.
3. Rubber dam punch with the function of adjusting the hole diameter.
4. The set of high elasticity clamps winged and wingless with excellent shape memory, which help to stabilize rubber dam on the tooth.

INTENDED USE

This medical device is intended for use together with rubber dam during dental treatment to isolate the working area. It enables to place rubber dam in oral cavity in right position.

- Frames stabilize the stretched RUBBER-DAM in order to protection of soft tissues. Rubber dam should be situated and stretched on the frame the way to avoid any unevennesses.
- Forceps with the opening degree resistance are used for putting on and taking off clamps. Use the forceps to stretch the clamp so that its wings become tense.
- Rubber-Dam punch is used for cutting holes in the RUBBER-DAM. Adjust the hole diameter to fit the tooth the clamp is being fitted. The punch has 6 holes of different diameter. To isolate anteriors or premolars use the holes of small diameter. Bigger holes should be used to isolate molar teeth. The biggest hole is for the tooth with fixed clamp.
- Clamps are used to stabilize the RUBBER-DAM inside the oral cavity.

DESCRIPTION OF PARTICULAR CLAMPS	CLAMP NUMBER	CLAMP SUITABLE FOR TOOTH NUMBER
Flat jaws for small lower molars.	3 	36, 37, 38, 46, 47, 48
Festooned wings for large upper molars. Wings fit the gum curvature.	5 	16, 17, 18, 26, 27, 28
Festooned jaws for molar's roots.	W8A 	16, 17, 18, 26, 27, 28, 36, 37, 38, 46, 47, 48
Serrated jaws for lower right and upper left molars.	12A 	26, 27, 28, 46, 47, 48
Serrated jaws for lower left and upper right molars.	13A 	16, 17, 18, 36, 37, 38
Small and large jaw for anteriors.	210 	11, 12, 13, 21, 22, 23, 31, 32, 33, 41, 42, 43
Flat jaws for anteriors.	212 	11, 12, 13, 21, 22, 23, 31, 32, 33, 41, 42, 43
Rounded jaws for premolars.	27 	14, 15, 24, 25, 34, 35, 44, 45
Rounded and flat jaws for premolars.	29 	14, 15, 24, 25, 34, 35, 44, 45
Small jaws and high bow for anteriors and premolars.	00 	13, 14, 15, 23, 24, 25, 33, 34, 35, 43, 44, 45

RUBBER-DAM set INSTRUCTION FOR USE

1) In case of root canal treatment only the tooth being treated should be isolated.





During treatment of cavities isolating the operating field should be ensured. During posterior teeth treatment, for example, the tooth situated behind the one being treated should be isolated (if it is possible) and the isolation should continue to the opposite middle incisor.

All the anteriors need to be isolated while treating, say from premolar to premolar to have the access and to allow the comparison with the opposite side what is important during cosmetic procedures.

In pedodontics, the teeth which are treated should be isolated. Usually it concerns two or three teeth.

2) To select particular points on RUBBER-DAM, in which the hole should be made, the RUBBER-DAM template might be used (pic.1). Place Rubber-Dam on top of the template, select a particular tooth to be isolated (pic.2, pic.3).

To make a perfect hole in the Rubber-Dam punch the hole then lift the dam over the pin (pic.4).

3) Place and stretch RUBBER-DAM sheet on to a metal frame (pic.5).

4) Adjust the clamp to the tooth using forceps.

In case of working with clamps without wings squeeze the clamp through the bottom of the hole and position it around the tooth using forceps.

When working with clamps with wings the RUBBER-DAM is stretched on wings, then all is positioned on the tooth by using forceps.




After the rubber-dam is positioned, push the RUBBER-DAM off the wings.

Clamp should be always located on the last exposed tooth.

The clamp's bow must be always placed toward patient's throat (pic.6).

5) To stabilize RUBBER-DAM place RUBBER-DAM cord in the interdental space between the exposed tooth and the one protected with dam (pic.7).

6) If the seal is required RUBBER-DAM liquid might be used (pic.8).

HANDLING DENTAL TOOLS MADE OF STAINLESS STEEL

Special processes have been undertaken during the manufacture of the dental instruments to ensure corrosion-resistance. In spite of this, the durability of these dental instruments depends on their proper use and correct maintenance.

Cleaning, rinsing and washing

The packaged instruments are non – sterile. Before sterilisation these new instruments should be washed thoroughly with warm water and using detergent. Before each use these dental instruments must be sterilised.

Sterilisation

All (reusable) dental instruments must be disinfected, then washed out to remove detergent and then dried before sterilization. Ultrasonic cleansing is recommended in particular as the most effective sterilization.

Recommended sterilisation methods:

Method

Steam autoclave
Dry air
Chemical vapours

Recommended temperature and duration

121°C, for 20 minutes
160°C, from 60 to 120 minutes
132°C, for appr. 20 minutes

PRECAUTIONS

Adhere to the instructions included and the additional guidelines for cleaning and sterilisation of the instruments.

STORAGE

Dental instruments should be stored in temperatures ranging between 5-30°C.

Before use protect the instrument against contamination or damage of the material structure (cracks, bends, fractures, peeling, etc.).

Failure to adhere to the above stated instructions decreases instrument durability and reduces instrument lifetime.

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