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RUBBER-DAM set

RUBBER DAM Instruments

BEFORE USE, READ THE FOLLOWING INSTRUCTION CAREFULLY
 TO BE USED BY DENTIST ONLY
 NON-STERILE - STERILIZE DIRECTLY BEFORE USE

SET INCLUDED:

1. Profiled rubber dam frames available in 3 sizes (S, M, L).
2. Rubber dam forceps for clamps placement with the opening degree resistance.
3. Rubber dam punch with the function of adjusting the hole diameter.
4. The set of high elasticity clamps winged and wingless with excellent shape memory (10 pcs.).
5. Base for clamps.

INTENDED USE

This medical device is intended for use together with rubber dam during dental treatment to isolate the working area. It enables to place rubber dam in oral cavity in right position.


- Frames stabilize the stretched RUBBER-DAM in order to protect the soft tissues. Rubber dam should be situated and stretched in a way avoiding any unevennesses.
- Forceps with the opening degree resistance are used for putting on and taking off clamps. Use the forceps to stretch the clamp so that its wings become tense.
- Rubber-Dam punch is used for cutting holes in the RUBBER-DAM. Adjust the hole diameter to fit the tooth the clamp is being fitted on. The punch has 6 holes of different diameter. To isolate anteriors or premolars use the holes of small diameter. Bigger holes should be used to isolate molar teeth. The biggest hole is for the tooth with fixed clamp.
- Clamps are used to stabilize the RUBBER-DAM inside the oral cavity.

DESCRIPTION OF PARTICULAR CLAMPS	CLAMP NUMBER	CLAMP SUITABLE FOR TOOTH NUMBER
Small and large jaw for anteriors.	210 	11, 12, 13, 21, 22, 23, 31, 32, 33, 41, 42, 43
Flat jaws for anteriors.	212 	11, 12, 13, 21, 22, 23, 31, 32, 33, 41, 42, 43
Small jaws and high bow for anteriors and premolars.	00 	13, 14, 15, 23, 24, 25, 33, 34, 35, 43, 44, 45
Rounded jaws for premolars.	27 	14, 15, 24, 25, 34, 35, 44, 45
Rounded and flat jaws for premolars.	29 	14, 15, 24, 25, 34, 35, 44, 45
Festooned jaws for molar's roots.	W8A 	16, 17, 18, 26, 27, 28, 36, 37, 38, 46, 47, 48
Flat jaws for small lower molars.	3 	36, 37, 38, 46, 47, 48
Festooned wings for large upper molars. Wings fit the gum curvature.	5 	16, 17, 18, 26, 27, 28
Serrated jaws for lower right and upper left molars.	12A 	26, 27, 28, 46, 47, 48
Serrated jaws for lower left and upper right molars.	13A 	16, 17, 18, 36, 37, 38

Additionally, there are clamps available for sale:

Wings, flat and large jaws for premolars.	2A 	14, 15, 24, 25, 34, 35, 44, 45
Festooned wings, flat jaws for large premolars.	208 	14, 15, 24, 25, 34, 35, 44, 45

RUBBER-DAM set INSTRUCTION FOR USE

HANDLING DENTAL TOOLS MADE OF STAINLESS STEEL

Special processes have been undertaken during the manufacture of the dental instruments to ensure corrosion-resistance. In spite of this, the durability of these dental instruments depends on their proper use and correct maintenance.

Cleaning, rinsing and washing

The packaged instruments are non-sterile. Before sterilisation these new instruments should be washed thoroughly with warm water and using detergent. Before each use these dental instruments must be sterilised.

Sterilisation

All (reusable) dental instruments must be disinfected, then washed out to remove detergent and then dried before sterilization. Ultrasonic cleansing is recommended in particular as the most effective sterilization. Instruments should be sterilized in accordance with national requirements.

Recommended sterilisation methods:

Method

Steam autoclave

Minimal parameters:

134°C, for 5 minutes

PRECAUTIONS

Adhere to the instructions included and the additional guidelines for cleaning and sterilisation of the instruments.

STORAGE

Dental instruments should be stored in temperatures ranging between 5-30°C.

Before use protect the instrument against contamination or damage of the material structure (cracks, bends, fractures, peeling, etc.).

Failure to adhere to the above stated instructions decreases instrument durability and reduces instrument lifetime.