

CERKAMED

MADE IN POLAND

MEDICAL COMPANY



Seize the impossible!

+ single-tool system
for every type of root canal

+ precise grip of the
broken instrument

+ thinnest working tip
on the market - only $\varnothing 0,3$ mm



+ the body made
of surgical steel

+ life-long warranty

About us

Polish

PRODUCTION



CERKAMED

We are **Polish producer** of dental materials



Products sales to over **100 countries!**




Cerkamed Team

btrPen

Broken
Tool
Remover

Dental tool for removing broken canal instruments.



Seize the impossible!

+ single-tool system
for every type of root canal

+ precise grip of the
broken instrument

+ thinnest working tip
on the market - only \varnothing 0,3 mm

+ the body made of surgical steel

+ life-long warranty



btr-pen.com

Only one on the market!

The tool combines new technology, fine quality and effective work. BTR Pen makes the removal of broken dental instruments a predictable, universal and repeatable procedure.

The usage of the ultrathin and highly-elastic working tip, together with ultimate-endurance nitinol loop, makes it easy to get to the broken canal instrument. The tool allows certain grip and safe removal of the broken canal instrument regardless of its position in the canal.

cerkamed.com

btrPen Broken Tool Remover

Included:



BTR Pen:

- ergonomic and comfortable handle equipped with grooves facilitating a secure hand grip
- made from high quality surgical stainless steel
- universal size
- designed for multiple use
- safe for autoclave sterilisation



Changeable tip with nitinol loop:

- the combination of a thin and highly elastic needle of Ø 0,3 mm with nitinol loop
- easiness in insertion into narrow and bent canals
- primary shape memory
- maximum endurance to stretching
- biocompatibility, corrosion resistance




Spring:

- full control of loop's tension and size
- optimized elasticity parameter
- guarantee of light and flawless moves
- freedom of fitting the spring toughness



Calibrator:

- quick and easy formation of the loop into desired size
- made from stainless steel
- designed for multiple use



NEW!

Endo Removal System

UPGRADED by K. Gończowski DDS, PhD, MSc


The ERS system is a universal set of instruments dedicated to the removal of broken parts of endodontic instruments from root canals under visual control by surgical microscope.

All parts of ERS are designed for **multiple use**



Handle

- one universal screwable handle matching all the types of working tips
- light construction and ergonomical design
- universal for both the right-handed and the left-handed dentists
- adapted to work in every position of operator's hand grip
- ergonomic trigger lever enables to regulate the strength of the grip on broken instrument



Working tip with **the movable wedge**

- available in 2 sizes (March 2019): 1.0 mm and 0.8 mm
- equipped with permanently joined luer-lock thread to perform effective disinfection and cleaning
- simple design resistant to user errors and mechanical damage
- microwindows with a diameter of 1 mm / 0.4 mm give perfect view under the microscope and serve as the additional mechanical retention!
- special design of the locking wedge effectively block inside the tube all types of the broken instruments including high taper rotary or reciprocal files
- smooth and extremely simple work thanks to the function of automatic reversing of the movable blocking wedge in the case of releasing the pressure on the trigger
- the tube can be inserted onto the broken instrument at an acute angle

Working tip with **ligature loop**

- available in 1 size of 0,7 mm
- the ligature wire has shape memory
- enables to grip the broken instrument both from the front and from the side of the closing cone
- strong spring which guarantees smooth and extremely simple work



Additional elements available for purchase

The starting set includes:

- handle
- working tip with the movable wedge of 0,8 mm

The full set includes:

- handle
- working tip with the movable wedge of 0,8 mm
- working tip with the movable wedge of 1,0 mm
- working tip with ligature loop of 0,7 mm



■ 017-008
Endodontic
micro lever



■ 015-195
Endodontic
micro probe
- bent

www.endo-rs.com



Dent Fresh ^{KIT} Smart Original

Specialised **mouthwash kit**

6 in 1

REFRESHES • CLEANS • PROTECTS • STRENGTHENS • SOOTHES • CARES

Perfect solution:

- ✓ In a trip
- ✓ For a business meeting
- ✓ Before a date
- ✓ During work
- ✓ After every meal

NEW in our offer!



be **SMART**

Dent Fresh Smart

it is the only one on the market, specialised disposable dose of mouthwash. Combination of a mouthwash, dental floss and interdental brush that guarantees complex care of the oral cavity and refreshed breath.



AVAILABLE PACKAGES: Dent Fresh Smart 12x15ml

Dent Fresh Smart KIT 10x15ml
+ tooth cleaners small
+ tooth cleaners big



Available packages:
- bottle of 200 g
- bottle of 400 g

! use with: ENDO-PREP CREAM



▲ CHLORAXID LIQUID for root canals rinsing
Sodium hypochlorite 2% or 5,25% in liquid



▼ CHLORAXID EXTRA LIQUID for root canals rinsing
Sodium hypochlorite 2% or 5,25% in liquid + surfactants



Available packages:
- bottle of 200 g
- bottle of 400 g

! use with: ENDO-PREP CREAM






▼ **CHLORAXID GEL** Preparation for root canals rinsing
Sodium hypochlorite 2% or 5,25% in gel



Available package:
- syringe 2 ml

▲ **CHLORAXID GEL EXTRA** Preparation for root canals rinsing
Sodium hypochlorite 2% or 5,25% in gel + surfactants





ECONOMIC package!



Available packages:
 - syringe of 2 ml, 5 ml, 10 ml
 - MEGA PACK - 3 x syringe 10 ml

! recommended to use with liquids from CHLORAXID line



▲ ENDO-PREP CREAM CREAM for root canals preparation EDTA 15% in cream, urea peroxide



Safety Adapter System™



Available packages:
 - bottle of 50 ml, 120 ml, 200 g

▲ ENDO-SOLUTION LIQUID for root canals widening EDTA 17% in liquid



Available packages:
- syringe of 2 ml, 5 ml, 10 ml

▲ ENDO-PREP GEL GEL for root canals preparation
EDTA 17% in gel



DLS™
DISPENSER LABELLING SYSTEM



Available package:
- bottle of 120 ml

ENDO-SOLUTION PREMIUM ▲ LIQUID for root canals widening
EDTA 17% in liquid + surfactants




Available packages:
- bottle of 200 g
- bottle of 400 g

use with: **HYDROCAL**



▲ **GLUCO-CHEX 2%** LIQUID for root canals rinsing
Chlorhexidine digluconate 2% in liquid



Available packages:
- bottle of 200 g
- bottle of 400 g

▲ **CITRIC ACID 40%** LIQUID for root canals rinsing



Available packages:
- bottle of 200 g

▲ **ISOPROPYL ALCOHOL**



Available packages:
- 100 pcs., 50 pcs., 20 pcs.

▲ ENDO-TOP Endo irrigation needles



Available for products:
CHLORAXID 2% EXTRA,
CHLORAXID 5,25% EXTRA,
CHLORAXID 2%,
CHLORAXID 5,25%,
GLUCO-CHEX 2%,
ENDO-SOLUTION,
ENDO-SOLUTION PREMIUM,
CITRIC ACID 40%



20 x + 20 x

▲ ENDO-SET LABELLED luer lock dispensers of 5 ml - 20 pcs + 20 pcs of endo irrigation needles ENDO-TOP

DLS™
DISPENSER LABELLING SYSTEM

Available for products:
CHLORAXID 2% EXTRA,
CHLORAXID 5,25% EXTRA,
CHLORAXID 2%,
CHLORAXID 5,25%,
GLUCO-CHEX 2%,
ENDO-SOLUTION,
ENDO-SOLUTION PREMIUM,
CITRIC ACID 40%



20 x

▲ ENDO-PACK LABELLED luer lock dispensers of 5 ml - 20 pcs



Available packages:
- syringe of 2,1 g
- MEGA PACK – 4 x syringe 2,1 g



Available packages:
- syringe of 2,1 g
- MEGA PACK – 4 x syringe 2,1 g



▲ **CALCIPAST** Calcium hydroxide paste
Calcium hydroxide

▲ **CALCIPAST+I** Calcium hydroxide in paste with iodoform
Calcium hydroxide, iodoform



Available package:
- bottle of 45 ml

▲ **CANAL CLEAN** LIQUID for root canals rinsing and drying
Isopropyl alcohol, acetone



Available packages:
- syringe of 2,1 g
- MEGA PACK – 4 x syringe 2,1 g

▲ **CALCIPAST FORTE** Calcium hydroxide in paste with iodoform and chlorophenol
Calcium hydroxide, iodoform, chlorophenol



perfect X-ray contrast





Available packages:
- syringe of 5 ml or 10 ml

▲ GLUCO-CHEX 2% GEL

Gel for root canals preparation
Chlorhexidine digluconate 2%



Available package:
- container of 10 g



▲ HYDROCAL Calcium hydroxide

Available package:
- container of 30 g



▲ IODOFORM *Iodoform pure p.a.*



The best in the world!



short setting time

plasticity



biocompatibility

strength

Available packages:

- MINI 3 x 0,14 g BIO MTA+ powder, 1 ml of BIO MTA+ liquid
- STANDARD 6 x 0,14 g BIO MTA+ powder, 1 ml of BIO MTA+
- MAXI 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml of BIO MTA+ liquid
- PRO 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml of BIO MTA+ liquid, MATRIX MTA+

BIO MTA+ ▲ Biomaterial for filling and rebuilding of root canals



Available packages:

- MINI 3 x 0,14 g MTA+ powder, 1 ml liquid
- STANDARD 6 x 0,14 g MTA+ powder, 1 ml liquid
- MAXI 10 x 0,14 g MTA+ powder, 1 ml liquid
- PRO 10 x 0,14 g MTA+ powder, 1 ml liquid, MATRIX MTA+

MTA+ ▲ Material for filling and rebuilding of root canals





- Possible variants:
- MTA+ Aplikator Premium – bent 0,8 mm
 - MTA+ Aplikator Premium – bent 1,2 mm
 - MTA+ Aplikator Premium – straight 0,8 mm
 - MTA+ Aplikator Premium – straight 1,2 mm

MTA+ APPLICATOR PREMIUM ▲

Reusable dental tools for handling and application of prepared measure of MTA+ material



- Possible variants:
- inner dimension of the piston 0,8 mm
 - inner dimension of the piston 1,2 mm

MTA+ APPLICATOR ▲

Reusable dental tools for handling and application of prepared measure of MTA+ material





DENTAL TOOLS FOR HANDLING AND APPLICATION OF MTA+ ▲

Reusable dental tools for handling and application of prepared measure of MTA+ material



▲ **MATRIX MTA+** Dosing block MTA+



★ SYNTEX ★

- A REAL KING AMONG ALL THE SEALERS



Available package:
 - syringe of 10 g of the preparation
 + mixing tips
 + mixing papers



▲ SYNTEX Root canals filling and sealing material

Color chart and sizes:

08	10								
15	45	90							
20	50	100							
25	55	110							
30	60	120	F1	F2					
35	70	130	F3	F4					
40	80	140	F5						

2%, 4%, 6%, ProTaper

HAND ROLLED



excellent quality

full range of sizes




useful round package

▲ GUTTA PERCHA POINTS



▲ GUTTA PERCHA PLUGGER Reusable dental tool for vertical condensation of gutta-percha



Available package:
bottle of 10 ml



ORANGE GUTTANE ▲ Oil for dissolving gutta-percha
Orange oil 100%



EUCALYPTOL ▲ Oil for dissolving gutta-percha
Eucalyptus oil 100%



Color chart and sizes:

08	10	2%, 4%, 6%, ProTaper	
15	45		
20	50	100	
25	55	110	
30	60	120	F1
35	70	130	F2
40	80	140	F3
			F4
			F5



ABSORBENT PAPER POINTS ▲

Available package:
bottle of 10 ml



METHYLENE CHLORIDE ▲ *Dichloromethane pure p.a.*




ENDO-ASPIRATOR ▲ for removal liquids from the root canal by aspirating



Available packages:
syringe 2 ml, 10 ml, 50 ml
MEGA PACK 3 x 10 ml
MONSTER PACK: 10 x 2 ml



Available packages:
syringe 2 ml, 10 ml, 50 ml
MONSTER PACK: 10 x 2 ml



Choose the one according to your preferences

BLUE ETCH ▲ Dental etchant
0-phosphoric acid 36%

▲ BLUE ETCH FLOW Flowable dental etchant
0-phosphoric acid 36%



Available packages:
YELLOW PORCELAIN ETCH syringe 2 ml
+ SILANE syringe 2 ml

SET: YELLOW PORCELAIN ETCH + SILANE ▲ Hydrofluoric acid 9,5%, methacrylic silane



Fractured porcelain



Fracture repaired with composite (Gradia)



Available packages:
syringe 2 ml

YELLOW PORCELAIN ETCH ▲ Porcelain etch Hydrofluoric acid 9,5%




Available packages:
syringe 2 ml

▲ SILANE Methacrylic silane



Available packages:
 - 2 ml syringe + applicators
 - Mega Pack 4 x 2 ml syringe + applicators



▲ INOX Gel promoting complete polymerisation of dental composites
 propanetriol (glycerol), Ph Eur purified water, gel base



Available packages:
 - PLS mini: 3 x 1 ml of preparation + set of applicators
 - PLS standard: 6 x 1 ml of preparation + set of applicators
 - PLS maxi: 10 x 1 ml of preparation + set of applicators



SPECIAL Gift



Available package:
bottle of 10 g

▲ **ALUSTAT** LIQUID for staunching the bleeding
Aluminium chloride 25%



Available packages:
syringe 5 ml, 10 ml
MEGA PACK 3 x 10 ml

▲ **ALUSTAT GEL** GEL for staunching bleeding
Aluminium chloride 25%

Available packages:
syringe 0,8 g
MEGA PACK 4 x 0,8 g



▲ **ALUSTAT FOAM** FOAM for staunching bleeding
Aluminium chloride 20%



Available package:
dropping bottle 10 ml
+ applicators

▲ **TCA DENT** Preparation for local cauterisation the mucous membrane
Trichloroacetic acid 35%

000.0
Smallest
size
retraction cord
Available on the
market

6 SIZES available:

- 000.0**
- 000**
- 00**
- 0**
- 1**
- 2**



**unique smooth
weave**

BEST-CORD NANO
Knitted retraction cords




*length matters
- 254 cm*



3 TYPES



RETRACTION CORD PACKERS Tools for packing retraction cords in the gingival pockets



Available packages:
syringe 2 ml
bottle of 50 ml
MONSTER pack: 10 x 2 ml

▲ **GLUCOSITE** LIQUID and GEL for gingival pockets
Chlorhexidine digluconate, hydrogen peroxide



Available packages:
dropping bottle 10 ml, 20 ml

▲ **EUGENOL** Eugenol oil *Eugenol oil 100%*




Available package:
syringe 2 ml

CANAL DETECTOR ▲ Root canals orifices indicator
Methylene blue



Available package:
syringe 2 ml

▲ **RED DETECTOR** Caries detector *Rhodamine*



**Bonding
time 1-7 min.**



Available package:
container of 50 g

▲ **ZINC OXIDE**
Powder for making pastes

▲ **ZINC OXIDE
FAST**
Zinc oxide - fast setting formula



Two variants available:

- **mint**
- **original**

Available packages:

- Start Pack concentrate in syringe 50 ml
- + bottle
- + dispenser
- concentrate in syringe 50 ml



▲ DENT FRESH Specialised Mouthwash



Aroma:

- **mint**
- **orange**

Available package:
- spray 200 ml


▲ PULP SPRAY SPRAY for pulp vitality testing

ENDO-STATION GO

ENDO-STATION GO

Line of endodontic organizers


ENDO-STATION GO EASY 2




ENDO-STATION GO EASY 3



ENDO-STATION GO GUTTA BOX 2X6



ENDO-STATION GO GUTTA BOX 2X3



ENDO-STATION GO GUTTA BOX 12



ENDO-STATION GO available in colors:*

- green ● blue ● yellow ● orange ●
- orange neon ● transparent ○

* ENDO-STATION GO GUTTA BOX available only transparent

ENDO-STATION GO MINI

ENDO-STATION GO MINI

New!



ENDO-STATION GO MINI available in colors:

- green ● blue ● transparent ○

Excellent workplace organization

ENDO-STATION GO MINI

Advantages:

- ✓ modern
- ✓ universal use
- ✓ in three colors



- transparent
- green
- blue

1



2



3



4



5



ADDITIONAL ACCESSORIES

Available packages:

1. ELASTIC APPLICATORS 10 pcs or 20 pcs
2. GEL APPLICATORS 20 pcs or 50 pcs
3. BRUSH APPLICATORS 10 pcs or 20 pcs
4. APPLICATORS FOR ALUSTAT FOAM 10 pcs
5. APPLICATORS FOR SYNTEX (mixing tips) 5 pcs
6. PAPER FOR PRODUCT MIXING 30 pcs
7. APPLICATION PIPETTES 5 pcs

6



7



28 **ADDITIONAL ACCESSORIES**

Available variants:

- plain BLUE: medium 0,18 mm or heavy 0,22 mm
- mint GREEN: medium 0,18 mm or heavy 0,22 mm
- mint BLACK: medium 0,18 mm

Available package:
36 pcs.



▲ RUBBER DAM Dental RUBBER-DAM



▲ RUBBER DAM TEMPLATE

Available packages:
syringe of 2ml
MEGA PACK 4 x 1,2ml



▲ RUBBER-DAM LIQUID Light-Cured RUBBER-DAM liquid



Available packages:
20 pcs. and 36 pcs.

▲ RUBBER-DAM NAPKINS



3 SIZES:
fine, medium, thick

Package containing 1,2 m of rubber dam cord

▲ RUBBER-DAM CORD Dental dam stabilizing cord

Perfect picture!

75
pcs.
Anti-reflex

- universal
- single use
- flexible

▲ DENTAL CONTRASTOR



- Set available:
- rubber dam punch
 - clamps forceps
 - rubber dam frames (3 size)
 - the set of clamps
 - base for clamps



New!

each item of Rubber dam Set available

separately!

12 types clamps

Also available in satin and black version




- Flat jaws for small lower molars (No 3).
- Festooned wings for large upper molars.
- Wings fit the gum curvature (No 5).
- Festooned jaws for molar's roots (No W8A).
- Serrated jaws for lower right and upper left molars (No 12A).
- Serrated jaws for lower left and upper right molars (No 13A).
- Small and large jaw for anteriors (No 210).
- Flat jaws for anteriors (No 212).
- Rounded jaws for premolars (No 27).
- Rounded and flat jaws for premolars (No 29).
- Small jaws and high bow for anteriors and premolars (No 00).
- Wings, flat and large jaws for premolars (No 2A).
- Festooned wings, flat jaws for large premolars (No 208).



▲ RUBBER DAM SET Tools for rubber dam

New!

- Set available:
- rubber dam punch
 - clamps forceps
 - rubber dam frame (M size)
 - clamps No 12A, 13A, 208, 210




ADDITIONALLY RECEIVE *for Free!*

- RUBBER-DAM (medium plain)
- RUBBER-DAM TEMPLATE PREMIUM with pen
- RUBBER-DAM LIQUID
- Light-cured rubber dam liquid
- RUBBER-DAM NAPKINS
- RUBBER-DAM CORD
- Dental dam stabilizing cord (medium)

▲ RUBBER-DAM SET MINI Start Pack Tools for rubber dam

Perfect protection!



Packing includes
20 pcs of disposable pads

ONE SIZE



▲ **PROTECTIVE PADS FOR A PATIENT**



Available colours:

- green ●
- blue ●
- pink ●
- beige ●
- white ●
- transparent ●



PROTECTIVE SHIELDS FOR A PATIENTS ▲

- ✓ Plaque removal
- ✓ Gum massage
- ✓ Soft and flexible
- ✓ Reusable



New!

1 pcs.



12 pcs.

Now also available in single packaging

Available colours:

green 

blue 

pink 

beige 

white 

transparent 

▶ COMFORT
with set of foils

▶ COMFORT LIGHT
with set of foils

▲ PROTECTIVE VISORS

▶ PROTECTIVE VISORS
ATTACHED TO GLASSES

▲ ANTISTATIC FOILS SET of removable protective foils

Available packages:
- standard 5 pcs.
- relax 4 pcs.
- relax light 4 pcs.

- ★ anti-fog effect
- ★ cleans
- ★ protects

Available package:
- spray 30 ml



3 in 1

◀ FOG
BUSTER
for Protective Visor



▲ GLASSES FOR POLYMERIZATION





- orange
- blue
- pink
- yellow
- green
- purple



Available packages:
 - 6 x syringe 1 g (colour: orange, blue, pink, yellow, green, purple) + applicators
 - syringe 1 g (chosen colour) + applicators

colourful way to a healthy smile!

RAINBOW FLOW ▲ Flowable composite

colour palette for kids



3D
stickers





- Available packs:
- FLUOR DEFENDER mini:
3 x 1 ml of preparation + set of applicators
 - FLUOR DEFENDER standard:
6 x 1 ml of preparation + set of applicators
 - FLUOR DEFENDER maxi:
10 x 1 ml of preparation + set of applicators

FLUOR DEFENDER ▲ Protective enamel with fluoride for children
0,1% fluorosilane (1600 ppm F⁻), HEMA



hand puppet
for the youngest patients



Available packages:
1 pcs., 2 pcs. or 3 pcs.

leaflet
„ROOT CANAL
TREATMENT”



Available packages:
20 pcs., 50 pcs. or 100 pcs.





Certificate of Bravery

Available package:
MIX - 20 pcs.



Boxes for milk teeth

Available package:
5 pcs.





Bracelet



guessing games
for 'brave
little patients

My teeth brushing
calendar

New! →





Buy set of **24 boxes** of mouthwash and you will get

a **JUTE BAG AND A DEDICATED CUP**

GREEN DENT THERAPY

Herbal mouth therapy guaranteed by leading company in dentistry



SPECIAL GIFT!

Set contains 24 boxes of mouthwash – 3 pieces of each kind

the greatest expert of herbal medicine



Mateusz Senderski MSc Eng.
selector of GREEN DENT THERAPY herbs,
author of 8 handbooks about herbal medicine

Use herbs!



Relief

Supports healing after dental treatments

Characteristics:

Composition of 8 herbs soothing, calming irritations and breath refreshing. Recommended to be used alternatively after dental treatment and in case of bad breath (halitosis).

Composition:

mint leaf, alchemilla, chicory roots, chamomile flower, sage, hartsease, iceland lichen, hemp-nettle herb.



Hemostatic

Reducing gums bleeding

Characteristics:

Composition of 8 herbs with astringent and soothing properties. Recommended to aid treatment for reducing gums bleeding.

Composition:

oak bark, sage, chamomile flower, marigold blossom, capsella, clayweed leaf, yellow sweet clover herb, hemp-nettle herb.



Sept

Soothing mouth inflammation

Characteristics:

Composition of 14 herbs with soothing and astringent properties. Recommended to be used alternatively in case of trush, mouth inflammation and purulent gum inflammations.

Composition:

sage, hyssop, achillea, thyme, marjoram, common toad flox, narrowleaf plantain, alchemilla, marigold blossom, chamomile flower, eyebright, peppermint, oak bark, iceland lichen.



Scaling

Soothing ailments after scaling

Characteristics:

Composition of 8 herbs with soothing and calming properties. Recommended to be used in case of after scaling conditions and to prevents dental plaque build up.

Composition:

horse chestnut, oak bark, hernaria glabra, root of common chicory, chickweed, archangel, clayweed leaf, raspberry leaf.



Gum

Soothing inflamed gums

Characteristics:

Composition of 8 herbs with astringent and soothing properties. Recommended to aid treatment in case of inflamed gums.

Composition:

sage, cornflower, knapweed, clayweed leaf, rubus leaf, artcium root, althea, achillea, marigold blossom.



Fresh

Soothing halitosis

Characteristics:

Composition of 8 herbs soothing and calming irritations and breath refreshing. Recommended to be used alternatively in case of bad breath (halitosis) as also after dental treatments

Composition:

peppermint, lavender flower, anise fruit, sage, hyssop herb, coriander fruit, hartsease, clayweed leaf.



Prosthesis

For persons wearing dentures

Characteristics:

Composition of 10 herbs with soothing and refreshing properties. Recommended to be used alternatively in case of wearing dentures conditions.

Composition:

thyme, sage, alchemilla, herba herniariae, althea, calamus, smartweed, clayweed leaf, great mullein, tilla platyphyllos.



Aphtha Stop

Supports healing of aphtha

Characteristics:

Composition of 8 herbs with astringent and soothing properties. Recommended to aid treatment of aphtha.

Composition:

licorice, chamomile flower, lemon balm, alchemilla, flower of calendula, chicory roots, iceland lichen, sage.



NEW IN OUR OFFER!



Living Dental Therapy

The only therapy on the market that comprehensively supports the patient's immune system before and after dental procedures.

iving Dental Therapy BEFORE

Living Dental Therapy Before provides the body with essential vitamins and minerals, protects the patient from oxidative stress during the procedure, supports processes related to the proper functioning of the immune system and accelerates wound healing.



vitamin C

collagen production
healthy teeth and gums



calcium

healthy bones
and teeth

zinc

protects from oxidative
stress


iving Dental Therapy AFTER

Living Dental Therapy After – a 14-day therapy for use after a dental procedure, which with its perfectly matched formula, effectively promotes wound healing, ensures the proper condition of oral mucous membranes and proper functioning of the immune system.



gotu kola

supports wound
healing



calcium

biotin B7

affects the immune
system and skin

vitamin B12

proper functioning of the
immune system

vitamin B2

proper condition of
mucous membranes





NOTICE! ALL TRADEMARKS ARE THE PROPERTY OF CERKAMED COMPANY.

THIS INFORMATION IS ONLY COMMERCIAL INFORMATION AND IS NOT AN OFFER WITHIN THE MEANING OF THE CIVIL CODE.

NEW
APPLICATION

CERKAMED
SMART ORDER




APPLICATION **SMART ORDER** CERKAMED



Our application will make your purchase even easier!

- ✓ Shortens your way to our products.
- ✓ Convenient, works on PC, mobile, tablets
- ✓ Readable and intuitive
- ✓ Provides information about desired products on your market
- ✓ Makes shopping with us fast and easy
- ✓ All you need is internet connection

Your purchase is just few clicks away!



www.cerkamed.com



CONTACT WITH US

PPH CERKAMED
Wojciech Pawłowski

NIP: PL 8652048770
REGON: 431162810

ul. Kwiatkowskiego 1, 37-450 Szałowa Wola
woj. podkarpackie, Poland

Trade department: Phone: +48 15 840 08 09
Phone: +48 15 842 35 85, int. 809
Mobile: 601 278 010

e-mail: export@cerkamed.pl



www.cerkamed.com

