

PRODUCT CATALOG

Cerkamed

Made in Poland 

MEDICAL COMPANY

Limited Edition



New golden tip:

0.5

0.4

0.3



CERKAMED

MEDICAL COMPANY

We are
Polish producer
of dental materials



Stalowa
Wola





Products
sales to
over

100 countries!



Cerkamed Team



btrPen Broken
Tool
Remover

Dental tool for removing broken canal instruments

Golden
edition


NEW

**BTR PEN
GOLDEN EDITION**

Let the **luxury**
enter you dental
office

NEW

**CALIBRATOR
GOLDEN EDITION**



Seize the impossible!

Join the world's endodontic elite
with our a **limited edition** of BTR PEN



New tips with
new sizes now
available!




facebook.com/groups/btrpen



btr-pen.com



Dental tool for removing broken canal instruments.



BTR Pen Author,
Piotr Pawłowski DDS

Seize the impossible!

+ single-tool system
for every type of root canal

+ precise grip of the
broken instrument

+ thinnest working tip
on the market - only Ø0,2 mm

+ the body made of surgical steel

+ life-long warranty



btr-pen.com

Only one on the market!

The tool combines new technology, fine quality and effective work. BTR Pen makes the removal of broken dental instruments a predictable, universal and repeatable procedure.

The usage of the ultrathin and highly-elastic working tip, together with ultimate-endurance nitinol loop, makes it easy to get to the broken canal instrument. The tool allows certain grip and safe removal of the broken canal instrument regardless of its position in the canal.



facebook.com/groups/btrpen



btr-pen.com

Included:

BTR PEN:

- ergonomic and comfortable handle equipped with grooves facilitating a secure hand grip
- made from high quality surgical stainless steel
- universal size
- designed for multiple use
- safe for autoclave sterilisation

Spring:

- full control of loop's tension and size
- optimized elasticity parameter
- guarantee of light and flawless moves
- freedom of fitting the spring toughness

Calibrator:

- quick and easy formation of the loop into desired size
- made from stainless steel
- designed for multiple use

Changeable tip with nitinol loop:

- the combination of a thin and highly elastic needle of $\varnothing 0,2$ mm with nitinol loop
- easiness in insertion into narrow and bent canals
- primary shape memory
- maximum endurance to stretching
- biocompatibility, corrosion resistance

**Impossible
changes
into possible**

0,3 mm

0,4 mm

New tips with new sizes
now available

0,5 mm

**Changeable tip
with nitinol loop**

5 pcs:
- 0.2 mm (available soon)
- 0.3 mm
- 0.4 mm
- 0.5 mm

NEW!**btr**


Dental training

DIGITAL MICROSCOPE**ADSM201**

- MINI HDMI input
- MINI USB port
- MICRO SD card input
- light regulation
- LED lightning

**MICROSCOPE'S
INSTRUCTION FOR USE****BTR Dental Training
KIT****Ideal kit for training**the removal of broken
endodontic files from root canals

The set includes:



2

1

6

5

4

3

**BTR Pen WORK
METHODOLOGY**

work with BTR Pen step by step

ANTI-SLIP MAT

• stability during the work with the kit

TWO TRAINING TEETH

with pieces of broken endodontic files

HARD 3 deeply placed, blocked pieces of files
(requiring processing before removal)

EASY 3 loosely placed pieces of files



SILICONE TRAINING MASS

- transparent and backlit
- 8 pieces of broken endodontic files
- a space to place the training tooth

SILICONE MASS STAND

- allows a 3D rotation

Also available:

btr training gum

The set includes:

SILICONE TRAINING MASS

with pieces of broken endodontic files

SILICONE MASS STAND


- allows a 3D rotation

ANTI-SLIP MAT

TWO TRAINING TEETH

with pieces of broken endodontic files (easy/hard)

BTR Dental Training GUM



3

CHLORAXID

LIQUID for root canals rinsing

Sodium hypochlorite 2%
or 5,25% in liquid

Available packages:

CHLORAXID 2% 200 g
CHLORAXID 2% 400 g
CHLORAXID 5,25% 200 g
CHLORAXID 5,25% 400 g



! use with: ENDO-PREP CREAM

! use with: ENDO-PREP CREAM

**CHLORAXID EXTRA**

LIQUID for root canals rinsing

Sodium hypochlorite 2%
or 5,25% in liquid + surfactants

Available packages:

CHLORAXID 2% EXTRA 200 g
CHLORAXID 2% EXTRA 400 g
CHLORAXID 5,25% EXTRA 200 g
CHLORAXID 5,25% EXTRA 400 g



! use with: ENDO-PREP CREAM

! use with: ENDO-PREP CREAM



CHLORAXID GEL

Preparation for root canals rinsing

Sodium hypochlorite 2% or 5,25% in gel

Available packages:

CHLORAXID 2% GEL 2 ml
CHLORAXID 5,25% GEL 2 ml



CHLORAXID GEL EXTRA

Preparation for root canals rinsing

Sodium hypochlorite 2% or 5,25% in gel + surfactants

Available packages:

CHLORAXID 2% GEL EXTRA 2 ml
CHLORAXID 5,25% GEL EXTRA 2 ml

ENDO-PREP CREAM

CREAM for root canals preparation

EDTA 15% in cream, urea peroxide

Available packages:


ENDO-PREP CREAM 2 ml
ENDO-PREP CREAM 5 ml
ENDO-PREP CREAM 10 ml
ENDO-PREP CREAM MEGA PACK
3 x 10 ml



! recommended to use with liquids from CHLORAXID line



Endo-Prep Cream



ENDO-PREP GEL

GEL for root canals preparation

EDTA 17% in gel

Available packages:

ENDO-PREP GEL 2 ml
ENDO-PREP GEL 5 ml
ENDO-PREP GEL 10 ml

ENDO-SOLUTION

LIQUID for root canals widening

EDTA 17% in liquid

Available packages:
ENDO-SOLUTION 50 ml
ENDO-SOLUTION 120 ml
ENDO-SOLUTION 200 g



Endo-Solution



ENDO-SOLUTION PREMIUM

LIQUID for root canals widening

EDTA 17% in liquid + surfactants

Available packages:
ENDO-SOLUTION PREMIUM 120 ml



GLUCO-CHEX 2%

LIQUID for root canals rinsing
Chlorhexidine digluconate 2% in liquid

Available packages:

GLUCO-CHEX 2% 200 g
GLUCO-CHEX 2% 400 g



! use with: HYDROCAL




CITRIC ACID 40%

LIQUID for root canals rinsing

Available packages:

CITRIC ACID 40% 200 g
CITRIC ACID 40% 400 g



ISOPROPYL ALCOHOL

Available packages:


ISOPROPYL ALCOHOL 200 g



ENDO-SET


LABELLED luer lock
dispensers of 5 ml - 20 pcs
+ 20 pcs ENDO-PACK
of endo irrigation needles
ENDO-TOP

Available for products:
CHLORAXID 2% EXTRA,
CHLORAXID 5,25% EXTRA,
CHLORAXID 2%, CHLORAXID 5,25%,
GLUCO-CHEX 2%,
ENDO-SOLUTION,
ENDO-SOLUTION PREMIUM,
CITRIC ACID 40%



0,2 mm

VS



0,3 mm



NEW

MICRO-TOP APPLICATORS

0.2 applicators with side hole for liquids

Available packages:
20 pcs.



ENDO-TOP

Endo irrigation needles

Available packages:
ENDO-TOP 20 pcs.
ENDO-TOP 50 pcs.
ENDO-TOP 100 pcs.



ENDO-PACK

LABELLED luer lock dispensers of 5 ml - 20 pcs

Available for products:
CHLORAXID 2% EXTRA,
CHLORAXID 5.25% EXTRA,
CHLORAXID 2%, CHLORAXID 5.25%,
GLUCO-CHEX 2%,
ENDO-SOLUTION,
ENDO-SOLUTION PREMIUM,
CITRIC ACID 40%

20 x



+

20 x



DLSTM
DISPENSER LABELLING SYSTEM

20 x





CALCIPAST

Calcium hydroxide paste
Calcium hydroxide

Available packages:

CALCIPAST 2,1 g
CALCIPAST MEGA PACK 4 x 2,1 g



CALCIPAST+I

Calcium hydroxide in paste
with iodoform

Calcium hydroxide, iodoform

Available packages:

CALCIPAST + I 2,1 g
CALCIPAST + I MEGA PACK 4 x 2,1 g



CALCIPAST FORTE

Calcium hydroxide in paste
with iodoform and
chlorophenol

Calcium hydroxide, iodoform,
chlorophenol

Available packages:

CALCIPAST FORTE 2,1 g
CALCIPAST FORTE MEGA PACK 4 x 2,1 g



before



after

CALCIPAST



CANAL CLEAN

LIQUID for root canals rinsing
and drying

Isopropyl alcohol, acetone

Available packages:

CANAL CLEAN bottle of 45 ml





GLUCO-CHEX 2% GEL

Gel for root canals preparation
Chlorhexidine digluconate 2%

Available packages:

GLUCO-CHEX 2% GEL 5 ml
GLUCO-CHEX 2% GEL 10 ml



HYDROCAL

Calcium hydroxide

Available packages:
HYDROCAL 10 g



JODOFORM

Available packages:
container of 30 g





The best in the world!

• short setting time • plasticity • biocompatibility • strength



BIO MTA+

Biomaterial for filling and rebuilding of root canals

Available packages:

MINI 3 x 0,14 g BIO MTA+ powder, 1 ml liquid
STANDARD 6 x 0,14 g BIO MTA+ powder, 1 ml liquid
MAXI 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml liquid
PRO 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml liquid, MATRIX MTA+



BIO MTA+ mini

Contains:
- 3 x 0,14 g BIO MTA+ powder, 1 ml BIO MTA+ liquid

BIO MTA+ standard

Contains:
- 6 x 0,14 g BIO MTA+ powder, 1 ml BIO MTA+ liquid

BIO MTA+ maxi

Contains:
- 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml BIO MTA+ liquid

BIO MTA+ pro

Contains:
- 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml BIO MTA+ liquid, MATRIX MTA+

MTA+ mini

Contains:
- 3 x 0,14 g MTA+ powder, 1 ml liquid

MTA+ standard

Contains:
- 6 x 0,14 g MTA+ powder, 1 ml liquid

MTA+ maxi

Contains:
- 10 x 0,14 g MTA+ powder, 2 x 1 ml liquid

MTA+ pro

Contains:
- 10 x 0,14 g MTA+ powder, 2 x 1 ml liquid, MATRIX MTA+



MTA+

Material for filling and rebuilding of root canals

Available packages:

MINI 3 x 0,14 g MTA+ powder, 1 ml liquid
STANDARD 6 x 0,14 g MTA+ powder, 1 ml liquid
MAXI 10 x 0,14 g MTA+ powder, 2x 1 ml liquid
PRO 10 x 0,14 g MTA+ powder, 2x 1 ml liquid, MATRIX MTA+



MTA+ APPLICATOR

Reusable dental tools for handling and application of prepared measure of MTA+ material

Possible variants:
APPLICATOR FOR MTA+
0,8 mm
APPLICATOR FOR MTA+
1,2 mm





DENTAL TOOLS FOR HANDLING AND APPLICATION OF MTA+

Reusable dental tools for handling and application of prepared measure of MTA+ material



MATRIX MTA+



★ SYNTEX ★

- A REAL KING AMONG ALL THE SEALERS

SYNTEX

Root canals filling and sealing material

Available packages:

syringe of 10 g of the preparation
+mixing tips
+mixing papers



GUTTA PERCHA POINTS

Available packages:

GUTTA PERCHA 2% 120 pcs.
GUTTA PERCHA 2% 120 pcs. marked
GUTTA PERCHA 4%, 6% 60 pcs.
GUTTA PERCHA 4%, 6% 60 pcs. marked
GUTTA PERCHA Pro Taper 60 pcs.
GUTTA PERCHA Pro Taper 60 pcs. marked
GUTTA PERCHA 2% ROUND PACKAGE 240 pcs.
GUTTA PERCHA 2% ROUND PACKAGE 240 pcs. marked
GUTTA PERCHA 4%, 6% ROUND PACKAGE 240 pcs.
GUTTA PERCHA 4%, 6% ROUND PACKAGE 240 pcs. marked



HAND ROLLED

08	10	
15	45	90
20	50	100
25	55	110
30	60	120
35	70	130
40	80	140
		F1
		F2
		F3
		F4
		F5


Color chart and sizes:

GUTTA PERCHA PLUGGER

Reusable dental tool for vertical condensation of gutta-percha

Available packages:

GUTTA PERCHA PLUGGER 1 pc.
GUTTA PERCHA PLUGGER 4 pcs.





EUCALYPTOL

Oil for dissolving gutta-percha
Eucalyptus oil 100%

Available packages:
EUCALYPTOL 10 ml
pipettes for application 5 pcs.



ORANGE GUTTANE

Oil for dissolving gutta-percha
Orange oil 100%

Available packages:
ORANGE GUTTANE 10 ml
pipettes for application 5 pcs.



METHYLENE CHLORIDE

Dichloromethane pure p.a.

Available packages:
METHYLENE CHLORIDE 10 ml



ABSORBENT PAPER POINTS

Available packages:
ABSORBENT PAPER POINTS 2% 200 pcs.
ABSORBENT PAPER POINTS 2% 200 pcs. marked
ABSORBENT PAPER POINTS 4%, 6% 100 pcs.
ABSORBENT PAPER POINTS 4%, 6% 100 pcs. marked
ABSORBENT PAPER POINTS ProTaper 100 pcs.
ABSORBENT PAPER POINTS ProTaper 100 pcs. marked
ABSORBENT PAPER POINTS 2%
Round package 400 pcs.
ABSORBENT PAPER POINTS 2%
Round package 400 pcs. marked
ABSORBENT PAPER POINTS 4%, 6%
Round package 240 pcs.
ABSORBENT PAPER POINTS 4%, 6%
Round package 240 pcs. marked 240 pcs.

ENDO-ASPIRATOR

for removal liquids from
the root canal by aspirating



ENDO aspirator^{PRO}

HIT  Warsaw
Dental Medica
Show

NEW

ENDO-ASPIRATOR PRO

for removal of liquids from root canals

Possible variants:

ENDO-ASPIRATOR PRO
ENDO-ASPIRATOR PRO Golden Version

PACKAGE INCLUDES:

HANDLE - 1 pc
ELASTIC APPLICATOR - 10 pcs
CONNECTOR - 3 pcs
LUER LOCK HOSE - 5 pcs
MAGNET - 5 pcs

equipped with **elastic applicator** allowing access to narrow root canals

magnetic fitting to the unit

safe for **autoclave sterilization**

German stainless steel

reusable

suitable for **every type of suction device**

World's only solution!



ADDITIONAL ACCESSORIES

PACKAGE INCLUDES:

LUER LOCK HOSE - 5 pcs
CONNECTOR - 1 pc




BLUE ETCH FLOW


1



2



3



BLUE ETCH

Dental etchant

O-phosphoric acid 36%

Available packages:

BLUE ETCH 2 ml
BLUE ETCH 10 ml
BLUE ETCH 50 ml
BLUE ETCH MEGA PACK 3 x 10 ml
BLUE ETCH MONSTER PACK 10 x 2 ml



BLUE ETCH FLOW

Flowable dental etchant

O-phosphoric acid 36%

Available packages:

BLUE ETCH FLOW 2 ml
BLUE ETCH FLOW 10 ml
BLUE ETCH FLOW 50 ml
BLUE ETCH FLOW MONSTER PACK
10 x 2 ml

YELLOW PORCELAIN ETCH

Hydrofluoric acid 9.5%

Available packages:
Syringe 2 ml



SILANE

Methacrylic silane

Available packages:
Syringe 2 ml

SET: YELLOW PORCELAIN ETCH + SILANE

Hydrofluoric acid 9.5%, methacrylic silane

Available packages:
YELLOW PORCELAIN ETCH syringe 2 ml
+ SILAN syringe 2 ml



Fractured porcelain



Fracture repaired with composite (Gratia)

INOX

Gel promoting complete polymerisation of dental composites


propanetriol (glycerol),
Ph Eur purified water, gel base

Available packages:

INOX 2 ml
INOX MEGA PACK 4 x 2 ml



STOP
Inhibition



PROTECT LIGHT SEAL

Dentine protective enamel

1% potassium fluoride (3500 ppm F),
HEMA

Available packages:

MINI 3 x 1 ml PROTECT LIGHT SEAL
+ set of applicators
STANDARD 6 x 1 ml PROTECT LIGHT SEAL
+ set of applicators
MAXI 10 x 1 ml PROTECT LIGHT SEAL
+ set of applicators



Special Gift

PULP SPRAY

SPRAY for pulp vitality testing

Available packages:

PULP SPRAY 200 ml

Aroma:

- mint
- orange
- miętowy
- pomarańczowy





ALUSTAT FOAM

FOAM for staunching the bleeding

Aluminium chloride 20%

Available packages:

ALUSTAT FOAM 0,8 g

ALUSTAT FOAM MEGA PACK 4 x 0,8 g

ALUSTAT GEL

GEL for staunching the bleeding

Aluminium chloride 25%

Available packages:

ALUSTAT GEL 5 ml

ALUSTAT GEL 10 ml

ALUSTAT GEL MEGA PACK 3 x 10 ml



ALUSTAT

LIQUID for staunching the bleeding

Aluminium chloride 25%

Available packages:

ALUSTAT 10 g



TCA DENT

Preparation for local cauterisation the mucous membrane

Trichloroacetic acid 35%

Available packages:

dropping bottle 10 ml

+ applicators



000.0 smallest size
retraction cord available on the market

BEST-CORD NANO

Knitted retraction cords,
254cm lenght

unique smooth
weave pła się

6 SIZES:

000.0
000
00
0
1
2

3 TYPES

D=0,22 mm

D=0,20 mm

D=0,28 mm

1,5 mm

Ø1,5 mm

Ø3 mm

1

2

3

RETRACTION CORD PACKERS

Tools for packing retraction cords in the gingival pockets

Possible variants:

RETRACTION CORD PACKERS

1. UNIVERSAL
2. ROUND 1,5 mm
3. ROUND 3 mm

4. UNIVERSAL - GOLDEN VERSION
5. ROUND 1,5 mm - GOLDEN VERSION
6. ROUND 3 mm - GOLDEN VERSION



**GOLDEN
VERSION**

NEW





GLUCOSITE GEL

GEL for gingival pockets

Chlorhexidine digluconate,
hydrogen peroxide

Available packages:

GLUCOSITE GEL 2 ml
GLUCOSITE GEL 50 ml
GLUCOSITE GEL MONSTER PACK
10 x 2 ml



GLUCOSITE

LIQUID for gingival pockets

Chlorhexidine digluconate,
hydrogen peroxide

Available packages:

GLUCOSITE 2 ml
GLUCOSITE 50 ml
GLUCOSITE MONSTER PACK 10 x 2 ml




CANAL DETECTOR

Root canals orifices indicator

Methylene blue

Available packages:

CANAL DETECTOR 2 ml



RED DETECTOR

Caries detector

Rhodamine

Available packages:

RED DETECTOR 2 ml



ZINC OXIDE FAST

Zinc oxide - fast setting formula

Available packages:

ZINC OXIDE FAST 50 g



**Bonding time
1-7 min.**

ZINC OXIDE

Powder for making pastes

Available packages:

ZINC OXIDE 50 g

EUGENOL

Eugenol oil

Eugenol oil 100%

Available packages:

EUGENOL 10 ml
EUGENOL 20 ml



each item of
Rubber dam Set
available

separately!

12 types of clamps

Also available in satin and black version

RUBBER-DAM SET

Tools for rubber dam

Set available:


RUBBER-DAM PUNCH

CLAMPS FORCEPS

RUBBER-DAM FRAMES (3 sizes)

THE SET OF CLAMPS

BASE FOR CLAMPS



Flat jaws for small lower molars (No 3).
Festooned wings for large upper molars.
Wings fit the gum curvature (No 5).
Festooned jaws for molar's roots (No WBA).
Serrated jaws for lower right and upper left molars (No 12A).
Serrated jaws for lower left and upper right molars (No 13A).
Small and large jaw for anteriors (No 210).
Flat jaws for anteriors (No 212).
Rounded jaws for premolars (No 27).
Rounded and flat jaws for premolars (No 29).
Small jaws and high bow for anteriors and premolars (No 00).
Wings, flat and large jaws for premolars (No 2A).
Festooned wings, flat jaws for large premolars (No 208).



RUBBER-DAM SET MINI

Tools for rubber dam
Start Pack

Set available:

RUBBER DAM PUNCH

CLAMPS FORCEPS

RUBBER DAM FRAME (M SIZE)

CLAMPS NO 12A, 13A, 208, 210

ADDITIONALLY RECEIVE

for free!

- RUBBER-DAM (medium plain)
- RUBBER-DAM TEMPLATE PREMIUM with pen
- RUBBER-DAM LIQUID
Light-cured rubber dam liquid
- RUBBER-DAM NAPKINS
- RUBBER-DAM CORD
Dental dam stabilizing cord (medium)



NEW QUALITY!

NEW

RUBBER-DAM SET GOLDEN VERSION

Tools for rubber dam

Available:

RUBBER-DAM SET Golden Version

RUBBER-DAM clamp Golden Version

GOLDEN
VERSION

3	12A	210	00
5	13A	27	2A
W8A	210	29	208





RUBBER-DAM LIQUID

Light-Cured
RUBBER-DAM liquid

Available packages:

RUBBER-DAM LIQUID 1,2 ml
RUBBER-DAM LIQUID
MEGA PACK 4 x 1,2 ml

RUBBER-DAM

Dental RUBBER-DAM

Available package:
36 pcs.

Available variants:

- plain BLUE: medium 0,18 mm
or heavy 0,22 mm
- mint GREEN: medium 0,18 mm
or heavy 0,22 mm
- mint BLACK: medium 0,18 mm



RUBBER-DAM NAPKINS

Available packages:


RUBBER-DAM NAPKINS 36 pcs.

RUBBER-DAM TEMPLATE

Available variants:

RUBBER-DAM TEMPLATE PREMIUM
RUBBER-DAM TEMPLATE Silver Version
RUBBER-DAM TEMPLATE Golden Version

marker *Gratis*



NEW


RUBBER-DAM CORD

Dental dam stabilizing cord

Package containing
1,2m of rubber dam cord

Available packages:

RUBBER-DAM CORD fine
RUBBER-DAM CORD medium
RUBBER-DAM CORD thick



SILICONE RINGS

for marking of instruments
in stomatology

Available packages:

SILICONE RINGS 20 pcs.
SILICONE RINGS 100 pcs.



Available colours:

- black
- purple
- blue
- red
- green
- yellow
- white



new

elastic applicators

HIT **IDS**
2019


For:

liquids • creams • pastes • gels


APPLICATORS

27G, 25G, 22G, 21G


Available packages:
20 pcs.



ELASTIC APPLICATOR 27G
bent with two sideholes
length: 35 mm
diameter: 0,4 mm



ELASTIC APPLICATOR 27G
straight
length: 41 mm
diameter: 0,4 - 0,2 mm



ELASTIC APPLICATOR 27G
bent
length: 41 mm
diameter: 0,4 - 0,2 mm

ELASTIC APPLICATOR 25G
bent
length: 41 mm
diameter: 0,5 - 0,3 mm



ELASTIC APPLICATOR 25G
long bent
length: 56 mm
diameter: 0,5 - 0,3 mm



ELASTIC APPLICATOR 22G
curve
length: 32,3 mm
diameter: 0,7 - 0,5 mm



ELASTIC APPLICATOR 21G
long bent
length: 56 mm
diameter: 0,8 - 0,6 mm

ADDITIONAL ACCESSORIES

Available packages:

1. ELASTIC APPLICATORS
10 pcs.
ELASTIC APPLICATORS
20 pcs.

2. GEL APPLICATORS
20 pcs.
GEL APPLICATORS
50 pcs.





3. BRUSH APPLICATORS
10 pcs.
BRUSH APPLICATORS
20 pcs.

4. APPLICATORS FOR ALLUSTAT FOAM
10 pcs.

5. APPLICATORS FOR SYNTAX
5 pcs.

6. PAPER FOR PRODUCT MIXING
30 pcs.

7. APPLICATION PIPETTES
5 pcs.



Perfect protection!

PROTECTIVE PADS FOR A PATIENT

Packing includes:
20 pcs. OF DISPOSABLE PADS


ONE SIZE



PROTECTIVE SHIELDS FOR A PATIENTS

Available colours:

- green
- blue
- pink
- beige
- white
- transparent



SIMPLY TRAVEL SILICONE FLOSS

Silicone dental floss

- ✓ Plaque removal
- ✓ Gum massage
- ✓ Soft and flexible
- ✓ Reusable



Now also available
in single packaging

PROTECTIVE VISORS

Available colours:

- green
- blue
- pink
- beige
- white
- transparent

COMFORT LIGHT
with set of foils

COMFORT
with set of foils

PROTECTIVE VISORS
ATTACHED TO GLASSES

ANTISTATIC FOILS

SET of removable protective foils

Available packages:
4x RELAX (0,175 mm)
5x STANDARD (0,125 mm)
4x LIGHT RELAX (0,175 mm)

GLASSES FOR POLYMERIZATION



Perfect solution:

- ✓ In a trip
- ✓ For a business meeting
- ✓ Before a date
- ✓ During work
- ✓ After every meal

NEW

DENT FRESH SMART

Specialised Mouthwash KIT

Available packages:

DENT FRESH SMART ORIGINAL
12x15ml
DENT FRESH SMART KIT ORIGINAL
10x15ml



Dent Fresh Smart
12 x 15 ml



Dent Fresh Smart KIT
10x15 ml

+ tooth cleaners small
+ tooth cleaners big



6 in 1

DENT FRESH

Specialised Mouthwash

Available packages:

- Start Pack concentrate in syringe 50 ml
+ bottle
+ dispenser
- concentrate in syringe 50 ml

Two variants available:

- MINT
- ORIGINAL



be
SMART

**Dent Fresh Smart
Toothpaste**



is a gel toothpaste in compact and disposable packaging.
Small and handy, it will be very useful in every trip.
It provides daily care and protection of oral cavity.



SUGAR FREE



**ACTIVE FILM
TECHNOLOGY**



NEW

**DENT FRESH SMART
Toothpaste Original**

AVAILABLE PACKAGES:

Dent Fresh Smart
Toothpaste 14 x 2 ml

Dent Fresh Smart
Toothpaste 8 x 5 x 2 ml

RAINBOW FLOW

Flowable composite

Available packages:

- 6 x syringe 1 g (colour: orange, blue, pink, yellow, green, purple) + applicators
- syringe 1 g (chosen colour) + applicators



- orange
- blue
- pink
- yellow
- green
- purple



RAINBOW FLOW

Flowable composite



FLUOR DEFENDER

Protective enamel with fluoride for children

0.1% fluorosilan (1600 ppm F)
HEMA

Available packs:

- FLUOR DEFENDER MINI
3 x 1 ml of preparation + set of applicators
- FLUOR DEFENDER MAXI
10 x 1 ml of preparation + set of applicators


MARKETING MATERIALS FOR PATIENT

Available packages:


- 3D STICKERS 10 pcs., 20 pcs., 50 pcs.
- LEAFLET „ROOT CANAL TREATMENT“ DEDICATED FOR PATIENT 20 pcs., 50 pcs., 100 pcs.
- PUPPET 1 pcs., 2 pcs., 3 pcs.
- BOXES FOR MILK TEETH 5 pcs.
- DIPLOMA „BRAVE PATIENT“ mix 20 pcs.
- GUESSING GAMES FOR A BRAVE PATIENT 20 pcs.
- TEETH BRUSHING CALENDAR for kids
Set included: teeth brushing calendar, diploma, stickers, magnet
- PUZZLES FOR CHILDREN 36 ELEMENTS




ENDO-STATION GO
EASY 3




ENDO-STATION GO
GUTTA BOX 12




ENDO-STATION GO
GUTTA BOX 2X3



ENDO-STATION GO
EASY 2



ENDO-STATION GO
GUTTA BOX 2X6



ENDO-STATION GO

Line of endodontic organizers

Available in colors*:

- green
- blue
- yellow
- orange
- orange neon
- transparent

* ENDO-STATION GO GUTTA BOX available only transparent

ENDO-STATION GO MINI

Available in colors:


- green
- blue
- transparent

Perfect picture!



DENTAL CONTRASTOR

- universal
- single use
- flexible



Endo Removal System

UPGRADED by K. Gończowski DDS, PhD, MSc

The ERS system is a universal set of instruments dedicated to the removal of broken parts of endodontic instruments from root canals under visual control by surgical microscope.

All parts of ERS are designed for **multiple use**



Handle

- one universal screwable handle matching all the types of working tips
- light construction and ergonomical design
- universal for both the right-handed and the left-handed dentists
- adapted to work in every position of operator's hand grip
- ergonomic trigger lever enables to regulate the strength of the grip on broken instrument



NEW!

Ø 0,7



Ø 0,8



NEW!

Ø 0,55



Ø 1



Working tip with **the movable wedge**

- available in 4 sizes (outer diameter of the tube) : 1,0 mm ; 0,8 mm ; 0,7 mm ; 0,55 mm
- equipped with permanently joined luer-lock thread to perform effective disinfection and cleaning
- simple design resistant to user errors and mechanical damage
- microwindows with a diameter of 1 mm / 0,4 mm give perfect view under the microscope and serve as the additional mechanical retention! (**NOTE!!!** 0,55 mm - A microtip **WITHOUT** microwindows behind initial serration!)
- special design of the locking wedge effectively block inside the tube all types of the broken instruments including high taper rotary or reciprocal files
- smooth and extremely simple work thanks to the function of automatic reversing of the movable blocking wedge in the case of releasing the pressure on the trigger
- the tube can be inserted onto the broken instrument at an



017-008
Endodontic
micro lever



015-195
Endodontic
micro probe
- bent

www.endo-rs.com





NOTICE! ALL TRADEMARKS ARE THE PROPERTY OF CERKAMED COMPANY.




NEW
APPLICATION

CERKAMED
SMART ORDER



APPLICATION **SMART ORDER** CERKAMED

Responsive



Our application will make your purchase even easier!

- Shortens your way to our products.
- Convenient, works on PC, mobile, tablets
- Readable and intuitive
- Provides information about desired products on your market
- Makes shopping with us fast and easy
- All you need is internet connection

Your purchase is just few clicks away!



- FOLLOW US - #cerkamed



PPH CERKAMED

Wojciech Pawłowski



Address: ul. Kwiatkowskiego 1
37-450 Stalowa Wola
woj. podkarpacie, Poland



Phone: +48 15 840 08 09
Phone: +48 15 842 35 85, int. 809
Mobile: +48 601 278 010



E-mail: export@cerkamed.pl
WWW: www.cerkamed.pl

VAT: PL 8652048770
REGON: 431162810